

Product Datasheet

Human C-type Lectin Domain Family 2 Member D ELISA Kit (orb1213278)

Description	Human C-type Lectin Domain Family 2 Member D, CLEC2D ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CLEC2D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,849
Uniprot ID	Q9UHP7
Solubility (25°C)	4800ng/L
Sensitivity	9.59ng/L
Entrez	29121
Expiration Date	Please enquire.