

Product Datasheet

Human Complement Factor H-related Protein 5 ELISA Kit (orb1213282)

Description	Human Complement Factor H-related Protein 5, CFHR5 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	CFHR5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	64,419
Uniprot ID	Q9BXR6
Solubility (25°C)	96ng/mL
Sensitivity	0.14ng/mL
Entrez	81494
Expiration Date	Please enquire.