

Product Datasheet

Human Platelet Derived Growth Factor D ELISA Kit (orb1213286)

Description	Human Platelet Derived Growth Factor D, PDGFD ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	PDGFD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,848
Uniprot ID	Q9GZP0
Solubility (25°C)	1600ng/L
Sensitivity	3.14ng/L
Entrez	80310
Expiration Date	Please enquire.