

## Product Datasheet

# Human Cell Adhesion Molecule 1 ELISA Kit (orb1213318)

<b>Description</b>	Human Cell Adhesion Molecule 1, CADM1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-4500ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	CADM1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	48,509
<b>Uniprot ID</b>	<a href="#">Q9BY67</a>
<b>Solubility (25°C)</b>	4800ng/L
<b>Sensitivity</b>	9.53ng/L
<b>Entrez</b>	<a href="#">23705</a>
<b>Expiration Date</b>	Please enquire.