

Product Datasheet

Human C-C Motif Chemokine 7 ELISA Kit (orb1213324)

Description	Human C-C Motif Chemokine 7, CCL7 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CCL7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,200
Uniprot ID	P80098
Solubility (25°C)	2400ng/L
Sensitivity	5.36ng/L
Entrez	6354
Expiration Date	Please enquire.