

Product Datasheet

Human C-C Chemokine Receptor Type 9 ELISA Kit (orb1213328)

Description	Human C-C Chemokine Receptor Type 9, CCR9 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CCR9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,016
Uniprot ID	P51686
Solubility (25°C)	3200ng/L
Sensitivity	7.42ng/L
Entrez	10803
Expiration Date	Please enquire.