

## Product Datasheet

# Human An1-type Zinc Finger Protein 6 ELISA Kit (orb1213347)

<b>Description</b>	Human An1-type Zinc Finger Protein 6, ZFAND6 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15-3000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	ZFAND6
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	22,555
<b>Uniprot ID</b>	<b>Q6FIF0</b>
<b>Solubility (25°C)</b>	3200ng/L
<b>Sensitivity</b>	7.65ng/L
<b>Entrez</b>	<b>54469</b>
<b>Expiration Date</b>	Please enquire.