

Product Datasheet

Human 14-3-3 Protein Theta ELISA Kit (orb1213385)

Description	Human 14-3-3 Protein Theta, YWHAQ ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	YWHAQ
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,764
Uniprot ID	P27348
Solubility (25°C)	64ng/ml
Sensitivity	0.091ng/ml
Entrez	10971
Expiration Date	Please enquire.