

Product Datasheet

Human Neuropeptide S ELISA Kit (orb1213409)

Description	Human Neuropeptide S, NPS ELISA Kit
Reactivity	Human
Range	5-2000ng/L
Tested Applications	ELISA
Target	NPS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,103
Uniprot ID	P0C0P6
Solubility (25°C)	2400ng/L
Sensitivity	3.78ng/L
Entrez	594857
Expiration Date	Please enquire.