

Product Datasheet

Human Microtubule-associated Proteins 1A/1B Light Chain 3B ELISA Kit (orb1213414)

Description	Human Microtubule-associated Proteins 1A/1B Light Chain 3B, MAP1LC3B ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	MAP1LC3B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,688
Uniprot ID	Q9GZQ8
Solubility (25°C)	2400ng/L
Sensitivity	5.28ng/L
Entrez	81631
Expiration Date	Please enquire.