

Product Datasheet

Human Erythroferrone ELISA Kit (orb1213425)

Description	Human Erythroferrone, ERFE ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	ERFE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,279
Uniprot ID	Q4G0M1
Solubility (25°C)	24ng/ml
Sensitivity	0.014ng/ml
Entrez	151176
Expiration Date	Please enquire.