

Product Datasheet

Human Low Density Lipoprotein Receptor Related Protein 5 ELISA Kit (orb1213438)

Description	Human Low Density Lipoprotein Receptor Related Protein 5, LRP5 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	LRP5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	179,145
Uniprot ID	O75197
Solubility (25°C)	1600ng/L
Sensitivity	3.72ng/L
Entrez	4041
Expiration Date	Please enquire.