

Product Datasheet

Human 5-Hydroxytryptamine Receptor 7 ELISA Kit (orb1213453)

Description	Human 5-Hydroxytryptamine Receptor 7, 5-HTR7 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	HTR7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,555
Uniprot ID	P34969
Solubility (25°C)	4000ng/L
Sensitivity	10.69ng/L
Entrez	3363
Expiration Date	Please enquire.