

Product Datasheet

Human Synaptophysin ELISA Kit (orb1213486)

Description	Human Synaptophysin, SYP ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	SYP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,845
Uniprot ID	P08247
Solubility (25°C)	9600ng/L
Sensitivity	14.58ng/L
Entrez	6855
Expiration Date	Please enquire.