

Product Datasheet

Human Immunoglobulin Binding Protein 1 ELISA Kit (orb1213500)

Description	Human Immunoglobulin Binding Protein 1, IGFBP1 ELISA Kit
Reactivity	Human
Range	10-2000ng/ml
Tested Applications	ELISA
Target	IGBP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,222
Uniprot ID	P78318
Solubility (25°C)	2400ng/ml
Sensitivity	5.46ng/ml
Entrez	3476
Expiration Date	Please enquire.