

Product Datasheet

Human Interleukin-36 Gamma ELISA Kit (orb1213503)

Description	Human Interleukin-36 Gamma, IL1F9 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	IL36G
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,721
Uniprot ID	Q9NZH8
Solubility (25°C)	9600ng/L
Sensitivity	16.54ng/L
Entrez	56300
Expiration Date	Please enquire.