

Product Datasheet

Human Glutathione S-transferase Mu 1 ELISA Kit (orb1213534)

Description	Human Glutathione S-transferase Mu 1, GSTM1 ELISA Kit
Reactivity	Human
Range	0.1-40.mU/ml
Tested Applications	ELISA
Target	GSTM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,712
Uniprot ID	P09488
Solubility (25°C)	48mU/ml
Sensitivity	0.053mU/ml
Entrez	2944
Expiration Date	Please enquire.