

Product Datasheet

Human Complement Factor H-related Protein 2 ELISA Kit (orb1213543)

Description	Human Complement Factor H-related Protein 2, CFHR2 ELISA Kit
Reactivity	Human
Range	3-900ng/ml
Tested Applications	ELISA
Target	CFHR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,651
Uniprot ID	P36980
Solubility (25°C)	960ng/ml
Sensitivity	1.29ng/ml
Entrez	3080
Expiration Date	Please enquire.