

## Product Datasheet

### Human Defensin-6 ELISA Kit (orb1213597)

<b>Description</b>	Human Defensin-6, DEFA6 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-20ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	DEFA6
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	10,975
<b>Uniprot ID</b>	<a href="#">Q01524</a>
<b>Solubility (25°C)</b>	24ng/ml
<b>Sensitivity</b>	0.024ng/ml
<b>Entrez</b>	<a href="#">1671</a>
<b>Expiration Date</b>	Please enquire.