

Product Datasheet

Human Bcl-2-like Protein 1 ELISA Kit (orb1213598)

Description	Human Bcl-2-like Protein 1, BCLX ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	BCL2L1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,049
Uniprot ID	Q07817
Solubility (25°C)	96ng/mL
Sensitivity	0.13ng/mL
Entrez	598
Expiration Date	Please enquire.