

Product Datasheet

Human Glutathione S-transferase Omega-1 ELISA Kit (orb1213607)

Description	Human Glutathione S-transferase Omega-1, GSTO1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	GSTO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,566
Uniprot ID	P78417
Solubility (25°C)	4800ng/L
Sensitivity	9.16ng/L
Entrez	9446
Expiration Date	Please enquire.