

## Product Datasheet

# Human Interleukin 18 Receptor Accessory Protein ELISA Kit (orb1213609)

<b>Description</b>	Human Interleukin 18 Receptor Accessory Protein, IL-18RAP ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-20ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL18RAP
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	68,310
<b>Uniprot ID</b>	<a href="#">O95256</a>
<b>Solubility (25°C)</b>	24ng/ml
<b>Sensitivity</b>	0.022ng/ml
<b>Entrez</b>	<a href="#">8807</a>
<b>Expiration Date</b>	Please enquire.