

Product Datasheet

Human Sortilin ELISA Kit (orb1213616)

Description	Human Sortilin, SORT1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	SORT1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	92,068
Uniprot ID	Q99523
Solubility (25°C)	24ng/ml
Sensitivity	0.024ng/ml
Entrez	6272
Expiration Date	Please enquire.