

Product Datasheet

Human Trefoil Factor 2 ELISA Kit (orb1213633)

Description	Human Trefoil Factor 2, TFF2 ELISA Kit
Reactivity	Human
Range	20-4500ng/mL
Tested Applications	ELISA
Target	TFF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,284
Uniprot ID	Q03403
Solubility (25°C)	4800ng/mL
Sensitivity	10.28ng/mL
Entrez	7032
Expiration Date	Please enquire.