

Product Datasheet

Human Trefoil Factor 1 ELISA Kit (orb1213634)

Description	Human Trefoil Factor 1, TFF1 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	TFF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,150
Uniprot ID	P04155
Solubility (25°C)	160ng/ml
Sensitivity	0.22ng/ml
Entrez	7031
Expiration Date	Please enquire.