

## Product Datasheet

# Human Secreted Frizzled Related Protein 3 ELISA Kit (orb1213635)

<b>Description</b>	Human Secreted Frizzled Related Protein 3, SFRP3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-30ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	FRZB
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	36,254
<b>Uniprot ID</b>	<a href="#">Q92765</a>
<b>Solubility (25°C)</b>	40ng/ml
<b>Sensitivity</b>	0.021ng/ml
<b>Entrez</b>	<a href="#">2487</a>
<b>Expiration Date</b>	Please enquire.