

Product Datasheet

Human Insulin Like Growth Factor Binding Protein 7 ELISA Kit (orb1213651)

Description	Human Insulin Like Growth Factor Binding Protein 7, IGFBP7 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	IGFBP7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,130
Uniprot ID	Q16270
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	3490
Expiration Date	Please enquire.