

Product Datasheet

Human Endothelial Protein C Receptor ELISA Kit (orb1213656)

Description	Human Endothelial Protein C Receptor, EPCR ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	PROCR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,671
Uniprot ID	Q9UNN8
Solubility (25°C)	128ng/ml
Sensitivity	0.27ng/ml
Entrez	10544
Expiration Date	Please enquire.