

Product Datasheet

Human Bcl-2-like Protein 11 ELISA Kit (orb1213684)

| | |
|----------------------------|---|
| Description | Human Bcl-2-like Protein 11, BCL2L11 ELISA Kit |
| Reactivity | Human |
| Range | 15-3000ng/L |
| Tested Applications | ELISA |
| Target | BCL2L11 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 22,171 |
| Uniprot ID | O43521 |
| Solubility (25°C) | 3200ng/L |
| Sensitivity | 6.98ng/L |
| Entrez | 10018 |
| Expiration Date | Please enquire. |