

## Product Datasheet

### Human Protein S100-B ELISA Kit (orb1213691)

<b>Description</b>	Human Protein S100-B, S100B ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.5-300ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	S100B
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	10,713
<b>Uniprot ID</b>	<a href="#">P04271</a>
<b>Solubility (25°C)</b>	320ng/L
<b>Sensitivity</b>	0.22ng/L
<b>Entrez</b>	<a href="#">6285</a>
<b>Expiration Date</b>	Please enquire.