

## Product Datasheet

# Human Hepatoma Derived Growth Factor ELISA Kit (orb1213693)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Hepatoma Derived Growth Factor, HDGF ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 0.05-30ng/ml   |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | HDGF   |
| <b>Storage</b>             | 2-8°C  |
| <b>Note</b>                | For research use only                                |
| <b>Application notes</b>   | Sandwich   |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates               |
| <b>MW</b>                  | 26,788   |
| <b>Uniprot ID</b>          | <b>P51858</b>  |
| <b>Solubility (25°C)</b>   | 32ng/ml  |
| <b>Sensitivity</b>         | 0.031ng/ml   |
| <b>Entrez</b>              | <b>3068</b>  |
| <b>Expiration Date</b>     | Please enquire.                                      |