

Product Datasheet

Human Nicotinamide Phosphoribosyltransferase ELISA Kit (orb1213696)

Description	Human Nicotinamide Phosphoribosyltransferase, NAMPT ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	VISFATIN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,521
Uniprot ID	P43490
Solubility (25°C)	80ng/ml
Sensitivity	0.115ng/ml
Entrez	10135
Expiration Date	Please enquire.