

Product Datasheet

Human Actin ELISA Kit (orb1213713)

Description	Human Actin, Cytoplasmic 2, ACTG1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ACTG1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,793
Uniprot ID	P63261
Solubility (25°C)	3200ng/L
Sensitivity	7.92ng/L
Entrez	71
Expiration Date	Please enquire.