

Product Datasheet

Human Fatty Acid Synthase ELISA Kit (orb1213720)

Description	Human Fatty Acid Synthase, FAS ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	FASN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	273,427
Uniprot ID	P49327
Solubility (25°C)	128ng/ml
Sensitivity	0.21ng/ml
Entrez	2194
Expiration Date	Please enquire.