

Product Datasheet

Human Glutamyl Aminopeptidase ELISA Kit (orb1213746)

Description	Human Glutamyl Aminopeptidase, GLUAP ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	ENPEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	109,244
Uniprot ID	Q07075
Solubility (25°C)	64ng/ml
Sensitivity	0.09ng/ml
Entrez	2028
Expiration Date	Please enquire.