

Product Datasheet

Human Pro-Opiomelanocortin ELISA Kit (orb1213766)

| | |
|----------------------------|---|
| Description | Human Pro-Opiomelanocortin, POMC ELISA Kit |
| Reactivity | Human |
| Range | 30-9000ng/L |
| Tested Applications | ELISA |
| Target | POMC |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 29,424 |
| Uniprot ID | P01189 |
| Solubility (25°C) | 9600ng/L |
| Sensitivity | 14.84ng/L |
| Entrez | 5443 |
| Expiration Date | Please enquire. |