

Product Datasheet

Human Alpha-Crosslaps ELISA Kit (orb1213768)

Description	Human Alpha-Crosslaps, ACTX ELISA Kit
Reactivity	Human
Range	10-3000ng/L
Tested Applications	ELISA
Target	ENO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,169
Uniprot ID	P06733
Solubility (25°C)	3200ng/L
Sensitivity	5.45ng/L
Entrez	2023
Expiration Date	Please enquire.