

Product Datasheet

Human Phosphatidylinositol Binding Clathrin Assembly Protein ELISA Kit (orb1213784)

Description	Human Phosphatidylinositol Binding Clathrin Assembly Protein, PICALM ELISA Kit
Reactivity	Human
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	PICALM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,755
Uniprot ID	Q13492
Solubility (25°C)	80ng/ml
Sensitivity	0.092ng/ml
Entrez	8301
Expiration Date	Please enquire.