

Product Datasheet

Human Receptor Tyrosine-protein Kinase Erbb-2 ELISA Kit (orb1213785)

Description	Human Receptor Tyrosine-protein Kinase Erbb-2, ERBB2 ELISA Kit
Reactivity	Human
Range	0.1-38ng/ml
Tested Applications	ELISA
Target	ERBB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	137,910
Uniprot ID	P04626
Solubility (25°C)	40ng/ml
Sensitivity	0.043ng/ml
Entrez	2064
Expiration Date	Please enquire.