

## Product Datasheet

### Human Cd5 Antigen-like ELISA Kit (orb1213813)

<b>Description</b>	Human Cd5 Antigen-like, CD5L ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3-900ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	CD5L
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	38,088
<b>Uniprot ID</b>	<a href="#">O43866</a>
<b>Solubility (25°C)</b>	960ng/mL
<b>Sensitivity</b>	1.62ng/mL
<b>Entrez</b>	<a href="#">922</a>
<b>Expiration Date</b>	Please enquire.