

Product Datasheet

Human Myelin Associated Glycoprotein ELISA Kit (orb1213824)

Description	Human Myelin Associated Glycoprotein, MAG ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	MAG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,069
Uniprot ID	P20916
Solubility (25°C)	64ng/ml
Sensitivity	0.11ng/ml
Entrez	4099
Expiration Date	Please enquire.