

Product Datasheet

Human Transferrin Receptor Protein 2 ELISA Kit (orb1213847)

Description	Human Transferrin Receptor Protein 2, TFR2 ELISA Kit
Reactivity	Human
Range	5-1000ug/ml
Tested Applications	ELISA
Target	TFR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	88,755
Uniprot ID	Q9UP52
Solubility (25°C)	1280µg/ml
Sensitivity	2.29µg/ml
Entrez	7036
Expiration Date	Please enquire.