

Product Datasheet

Human Ficolin 1 ELISA Kit (orb1213877)

Description	Human Ficolin 1, FCN1 ELISA Kit
Reactivity	Human
Range	2-600ng/ml
Tested Applications	ELISA
Target	FCN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,078
Uniprot ID	O00602
Solubility (25°C)	800ng/ml
Sensitivity	0.92ng/ml
Entrez	2219
Expiration Date	Please enquire.