

## Product Datasheet

### Human Signal Peptide ELISA Kit (orb1213878)

<b>Description</b>	Human Signal Peptide, Cub and EGF-like Domain-containing Protein 1, SCUBE1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1-400ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	SCUBE1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	107,910
<b>Uniprot ID</b>	<a href="#">Q8IWY4</a>
<b>Solubility (25°C)</b>	480ng/ml
<b>Sensitivity</b>	0.55ng/ml
<b>Entrez</b>	<a href="#">80274</a>
<b>Expiration Date</b>	Please enquire.