

Product Datasheet

Human Neurofilament Heavy Polypeptide ELISA Kit (orb1213882)

Description	Human Neurofilament Heavy Polypeptide, NEFH ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	NEFH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	112,479
Uniprot ID	P12036
Solubility (25°C)	160ng/ml
Sensitivity	0.27ng/ml
Entrez	4744
Expiration Date	Please enquire.