

Product Datasheet

Human C-C Motif Chemokine 26 ELISA Kit (orb1213888)

Description	Human C-C Motif Chemokine 26, CCL26 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CCL26
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,648
Uniprot ID	Q9Y258
Solubility (25°C)	2400ng/L
Sensitivity	5.49ng/L
Entrez	10344
Expiration Date	Please enquire.