

Product Datasheet

Human Platelet-derived Growth Factor C ELISA Kit (orb1213900)

Description	Human Platelet-derived Growth Factor C, PDGF-C ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	PDGFC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,029
Uniprot ID	Q9NRA1
Solubility (25°C)	1280ng/L
Sensitivity	3.12ng/L
Entrez	56034
Expiration Date	Please enquire.