

Product Datasheet

Human Myelin-oligodendrocyte Glycoprotein ELISA Kit (orb1213901)

Description	Human Myelin-oligodendrocyte Glycoprotein, MOG ELISA Kit
Reactivity	Human
Range	10-2000ng/mL
Tested Applications	ELISA
Target	MOG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,193
Uniprot ID	Q16653
Solubility (25°C)	2400ng/mL
Sensitivity	4.84ng/mL
Entrez	4340
Expiration Date	Please enquire.