

Product Datasheet

Human Krueppel-like Factor 5 ELISA Kit (orb1213906)

Description	Human Krueppel-like Factor 5, KLF5 ELISA Kit
Reactivity	Human
Range	10-4000pmol/L
Tested Applications	ELISA
Target	KLF5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,792
Uniprot ID	Q13887
Solubility (25°C)	4800pmol/L
Sensitivity	5.21pmol/L
Entrez	688
Expiration Date	Please enquire.