

## Product Datasheet

# Human Stress Induced Phosphoprotein 1 ELISA Kit (orb1213920)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Stress Induced Phosphoprotein 1, STIP1 ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 5-1500ng/L   |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | STIP1  |
| <b>Storage</b>             | 2-8°C  |
| <b>Note</b>                | For research use only                                  |
| <b>Application notes</b>   | Sandwich   |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                 |
| <b>MW</b>                  | 62,639   |
| <b>Uniprot ID</b>          | <b>P31948</b>  |
| <b>Solubility (25°C)</b>   | 1600ng/L   |
| <b>Sensitivity</b>         | 2.42ng/L   |
| <b>Entrez</b>              | <b>10963</b>   |
| <b>Expiration Date</b>     | Please enquire.  |